CELLYARD[®] beads 200

Surface activated cell carrier beads
for high density cell culture

INSTRUCTION MANUAL


Fig.1 SEM photograph of CELLYARD beads


Fig.2 SEM photograph of the surface of nylon particle


Fig.3 SEM photograph of the surface of CELLYARD beads


Fig.4 The C6/36 cells adhering to the surface of CELLYARD beads

1. USE

The high-density mass culture of various attachment dependent cell .

2. MERIT

- 1) The CELLYARD beads (Fig.1) is composed of microspheres that hydroxyapatite $[\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2]$ fine particles cover the surface of the nylon particles of average particle size $200 \mu\text{m}$. Hydroxyapatite is one of calcium phosphate system compounds , has excellent biocompatibility , and has already been used as artificial bone , bone prosthetic material , and so on .

In comparison with surface of the nylon particle (Fig.2) and surface of the CELLYARD beads (Fig.3) , the structure of the nylon which exists on the surface of the nylon particle is not completely observed in the surface of the CELLYARD beads . This means that the particle surface of the CELLYARD beads is covered with hydroxyapatite nearly completely . Therefore , the CELLYARD beads has the capability to adsorb various protein , a cell , etc. like hydroxyapatite (Fig.4) .¹⁾

- 2) The particle size of the CELLYARD beads is suitable for adhesion and multiplication of cells .²⁾
- 3) The specific gravity of the CELLYARD beads ($1.01 \sim 1.05 \text{ g/cm}^3$) leads suspension in calm agitation and spontaneous sedimentation in standing condition , and is suitable for recovery of the CELLYARD beads or medium .
- 4) The CELLYARD beads do not swell in solution .

3. USAGE

The CELLYARD beads is suitable for high-density culture of various attachment dependent cell, but optimum condition, such as capacity of culture flask, form, and agitation speed, change with kind of cell. Therefore, this instruction manual explains the method for culture of the Vero cell (the African green monkey kidney origin cell) using spinner flask for 100 mL.

Fig.5 CELLYARD beads after sterilization

Fig.6 The suspension of CELLYARD beads is moved to centrifugal tube.

Fig.7 CELLYARD beads is suspended to the medium.

Fig.8 The recovery of the cell

1) STERILIZATION

50g the CELLYARD beads suspended to 500 mL PBS(-) are carried out the autoclave sterilization (121°C, 20 min). The CELLYARD beads is finely dispersed after sterilization, although it may be hard to disperse to PBS(-) at first. Sterilized suspension (100 g/l) is conserved at room temperature (Fig.5).

※ Please perform operation of following 2)–4) within the clean bench of aseptic condition.

2) WASHING

10 mL of suspension of the CELLYARD beads prepared in 1) is moved to another centrifugal tube which sterilized (Fig.6), and supernatant is discarded by centrifugation (1600 rpm, 5 min). You may carry out spontaneous sedimentation without using centrifuge. 10 mL of the MEM medium (FCS free) is added to the pellet of the CELLYARD beads, and is washed. Then, centrifugation (1600 rpm, 5 min) or spontaneous sedimentation is carried out and supernatant is discarded. Please repeat this operation 3 times.

10 mL of MEM medium (10% FCS) is newly added to the pellet of the CELLYARD beads in order to prepare the suspension of the CELLYARD beads.

※※ Please use a pipette with the large path of the hole at a tip.


Fig.9 Suspension of the CELLYARD beads and the cell are added to the spinner flask .

Fig.10 It is stood in the adhesion stage of the cell without agitation .

Fig.11 The medium is added to the final volume .


Fig.12 The CELLYARD beads in which the Vero cell adhered . (The culture start post-3 hour course)


Fig.13 The CELLYARD beads that the surface is perfectly covered by the Vero cell . (The culture start post-3 day course)

3) PREPARATION OF THE CELL

The cell is recovered at 20 mL of MEM medium (10% FCS) in order to prepare the suspension of the cell , so that the concentration of seeding cell may consist as 1.0×10^7 cells (Fig.8) .

4) SEEDING OF THE CELL

In sterilized spinner flask (for 100 mL) , 10 mL of suspension of the CELLYARD beads prepared in 2) and 20 mL suspension of the cell prepared in 3) are added , and the whole volume is done at about 30 mL , 1/3 of the final volume (Fig.9) .

5) CULTURE START

The spinner flask is stood in incubator (37°C, 5% CO₂) in the condition of 4) . It is stood for 3 hours (the adhesion stage of the cell) from the culture start without connected agitation (Fig.10) . In the time, please do 1 minute agitation at 37 rpm speed every 60 minute . In this time agitation speed , you do not mind at the speed in the case of the connected agitation in the vegetation stage of the cell . 3 hours passed, and MEM medium (10% FCS) is added in the clean bench to the final volume (100 mL) , and the connected agitation is started (Fig.11) .

The aspect of the cell which passed of culture start post-3 hour and for 3 days is shown in (Fig.12) and (Fig.13) .

4. ATTENTION IN THE HANDLING

- 1) By seeding cell, standing time in the adhesion stage must be adjusted in the about 3–6 hour. And, there is the case in which the connected agitation should be done by adding the medium of the final volume right after the seeding of the cell, and the agitation should not be completely done, while it is being stood in the adhesion stage. By adjusting optimum culture condition of seeding cell, please use the **CELLYARD beads**.
- 2) It is recommended that the medium used for washing of the **CELLYARD Beads** and culture are beforehand equilibrated in incubator at 37°C and 5%CO₂. So, the **CELLYARD beads** can be used more effectively.
- 3) Please do not dip the **CELLYARD beads** in acid solution and solution including chelating reagent and organic solvent such as acetone, xylene, and so on. The performance may deteriorate.
- 4) There is no problem on the performance, though the particle of black or brown may have been mixed in the product.
- 5) Please dispose of the **CELLYARD beads** as an incombustible refuse after the use.

5. REFERENCE

- 1) Yamamoto A., Nakajima T., Tominaga Y. and Ogawa T.: Porus hydroxyapatite ceramics, Tiss. Cult. Res. Commun. 18, 235–244, 1999
- 2) Shibata T.: The world of the bioreactor—the base and application for the practice person, Hario Laboratory Co., Ltd., 231, 1992
- 3) Sugo K., Kato M., Ishikawa T., Yamamoto A. and Ogawa T.: Hydroxyapatite microcarrier (1), Tiss. Cult. Res. Commun. 25, 105–111, 2006
- 4) Sugo K., Kato M., Ishikawa T., Yamamoto A. and Ogawa T.: Hydroxyapatite microcarrier (2), Tiss. Cult. Res. Commun. 25, 113–118, 2006
- 5) Sugo K. and Ogawa T.: Three-dimensional Culture of Rat Bone Marrow Cells Using Hydroxyapatite Microcarrier, Tiss. Cult. Res. Commun. 26, 125–131, 2007
- 6) Sugo K., Nakayama M., and Ogawa T.: Japanese Encephalitis Virus Production Using Hydroxyapatite Microcarrier, Medicine and Biology. 151 (7), 237–243, 2007

E200804

Distributor

COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618

CELLYARD™ beads

Surface activated cell carrier beads
for high density cell culture

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